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In vitro fertilization medium (Mouse-specific)

Product	Item No.	specification	Application
HTF	72002	2mL	Mouse fertilization fluid-Fresh Sperm
HTF-A	72102	2mL	Mouse fertilization fluid-cryorecovered sperm

Composition

Water for embryo transfer, BSA, Penicillin-Streptomycin, Phenol Red.etc

Uses and Principles

The reagent is mainly used for in vitro fertilization of fresh sperm and frozen-thawed sperm in mice, which generally needs to be used with c-TYH. With NOD strains,HTF can be used as sperm capacitation medium.

Storage conditions and expiration date

1. Storage: 4°C, Protect from light.
2. Expiration date: 6 months. If the reagent is not used after opening, please seal it in time and use it within 24 hours.

Mouse

Male mice in breeding colonies should be left alone for at least 3 days before use, and male mice in no- breeding colonies should be left alone for at least one week after trial mating before the experiment.

Select female of C57BL/6 (3.5-4 week’s old)strains ,BALB/CJ and NOD strains were recommended at 8 weeks of age.


Methods

1. IVF- (fertilization dish for fresh sperm)

1.1 The 200μL drops were used for fertilization,transfer cumulus masses from all three femals to a single fertilization drop. Incubate the fertilization dishes at 37°C.

1.2 Collect oocyte-cumulus complexes by tearing the oviduct ampulla under a stereos -copic microscope in a drop of 200μL HTF.

1.3 2-cells will be picked directly the next day.

HTF medium

2. IVF- (HTF could be use as sperm capacitation medium for NOD strains)

2.1 A 35-mm culture dish with HTF covered with oil was prepared and equilibrated for at least 30 min in the incubator.

2.2 Collect tail of epididymis from a male killed via cervical dislocation or anesthetized, and place them on a filter paper. Remove the fat and blood from the epididymides.

2.3 A 1 mL syringe was used to gently punctured 4-6 times in the tail of the epididymis, Squeeze tail of the epididymis by the hands and help the sperm outflow.

Transfer of sperm clumps using microforceps for placement in HTF for 60 min in the incubator .

2.4 After sperm capacitation, pick up 10 μ L of sperm suspension from the edge of the drop was added to the HTF dishes containing oocytes. Depending on the density and viability of the sperm, it is recommended to add about 8 μ L of sperm.


Usage - HTF-A

1. IVF-For preparation of cryorecovered sperm

1.1 Take 30.7 mg of A powder in an EP tube and add 1 mL of HTF solution, cover and mix gently to prepare "Solution A".

1.2 Pipette 10 μ L of "Solution A" into 990 μ L of HTF solution and mix well, this solution is called "Solution B" and will be used to make fertilized drops.

1.3 In a 35mm dish make a 90 μ L fertilization drop with "Solution B", then make three 50 μ L droplets and a 80 μ L drop with HTF, two 50 μ L droplets for washing and another 50 μ L droplet for culture for 2-cell culture, in addition, 80 μ L of HTF culture drop was used for the culture after washing. The specific functions are shown in the figure below.

1.4 The HTF dish was equilibrated in an incubator for 30 minutes before use. In general, one fertilization dish corresponds to one donor female mouse (about 30-40 eggs). Add 10 μ L of capacitated sperm to each fertilization drop(Solution B). After adding sperm for 4 hours, wash off granulosa cells and sperm, then culture in the 80 μ L drop of HTF.

Cautions

1. Avoid blood, fat and hair and tissue from entering capacitation drops during oocyte collection.
2. Solution B needs to be formulated when using, not to exceed 2 hours about time of use , otherwise it will affect the fertilization rate.
3. After sperm capacitation, it is recommended that pick up sperm suspension from the edge of the dropwas added to the HTF dishes.
4. HTF is packed in ampoule, the bottle is thin, please ensure safety when using it.