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**Sperm capacitation medium (Mouse-specific)**

| Product | Item No. | specification | Application        |
|---------|----------|---------------|--------------------|
| c-TYH   | 72021    | 2mL           | Sperm capacitation |

**Composition**

Water for embryo transfer, Penicillin-Streptomycin, Phenol Red,  $\beta$ -Cyclodextrin

**Storage conditions and expiration date**

1. Storage Temperature: 4°C, Protect from light.
2. Expiration date: 6 months. If the reagent is not used after opening, please seal it with sealing film in time and use it within 24 hours.

**Mouse**

12- to 24-week-old males, male mice in breeding colonies should be left alone for at least 3 days before use, and male mice in no- breeding colonies should be left alone for at least one week after trial mating before the experiment.

**Methods**

1. In vitro fertilization-fresh sperm
  - 1.1 A 35-mm culture dish with a 200  $\mu$ L drop of c-TYH covered with oil was prepared and equilibrated for at least 30 min in the incubator.
  - 1.2 Collect tail of epididymis from a male killed via cervical dislocation or anesthetized, and place them on a filter paper. Remove the fat and blood from the epididymides.
  - 1.3 A 1 mL syringe was used to gently punctured 4-6 times in the tail of the epididymis, Squeeze tail of the epididymis by the hands and help the "sperm clumps" outflow. Transfer of sperm clumps using microforceps for placement in the c-TYH for 30min in the incubator.
  - 1.4 After sperm capacitation, pick up sperm suspension from the edge of the drop was added to the HTF dishes containing oocytes.

Depending on the density and viability of the sperm, it is recommended to add about

5-8  $\mu\text{L}$  of sperm.

2. IVF-cryorecovered sperm

2.1 A 35-mm culture dish with a 90  $\mu\text{L}$  drop of c-TYH covered with oil was prepared and equilibrated for at least 30 min in the incubator.

2.2 Pick up the frozen straw from the liquid nitrogen tank, and place it in a 37 °C water bath for 10min.

2.3 Slowly aspirated sperm from straw and added it to the pre-equilibrated Sperm capacitation dish, capacitated in an incubator for 30min. pick up 10 $\mu\text{L}$  of sperm suspension from the edge of the drop was added to the HTF dishes containing oocytes.

**Cautions:**

1. Avoid blood, fat and hair and tissue from entering capacitation drops during Sperm capacitation medium.
2. Place the sperm capacitation dish in an incubator equilibrate for at least 30 min.
3. The angle of sperm aspiration needs to be inclined and gentle to avoid damaging the sperm tail.
4. Without tracing back the parental generation, if the sperm is low densities, the sperm of different strains mice can be collected into the same fertilization drop for improve the efficiency of IVF.