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Cryoprotective agent (Mouse-specific)

Product	Item No.	specification	Application
g-CPA	72022	1 mL	Sperm cryoprotectant solution in mice

Composition

Raffinose, skimmed milk powder, L-glutamine, etc.

Uses and Principles

This reagent is mainly used for sperm cryopreservation of mice, It is to be used with sperm straws, sealers, HTF, etc. The role of cryoprotective agent is to improve the motility and survival rate of sperm after frozen-thawed sperm. Reduces intracellular damage caused by change in osmolarity.

Storage conditions and expiration date

1. Storage: 4°C-25°C, Protect from light.
2. Expiration date: 6 months. If the reagent is not used after opening, please seal it with sealing film in time.

Mouse

Male mouse 3-4 month of age, male mice in breeding colonies should be left alone for at least 3 days before use, and male mice in no- breeding colonies should be left alone for at least one week after trial mating before the experiment.


Methods

1. Preparation of dishes: g-CPA was prepared and equilibrated for at least 30 min at 37°.A 35mm culture dish with a 120μL g-CPA drop (for washing the epididymis) , another 120μL drop of g-CPA covered with oil was placed on the heating plate for sperm suspension.
2. Harvesting vas deferens and cauda epididymis: Humanely euthanize male mice by cervical dislocation,Cut the vas deferens and gently pull out the vas deferens, removing most of the fat , hair and blood vessels.
3. Cut off the fat tissue just below the cauda epididymis and place the vas deferens and

cauda epididymis into the 35mm dish containing 120 μ L g-CPA (for washing the epididymis) .

4. Transfer cauda epididymis into another g-CPA drop (for release the sperm) ,use spring scissors to make several cuts for all cauda epididymis to release the sperm, Place the dish on 37 $^{\circ}$ C hot plate for 3 min. gently swirled every min for 20 sec to help the sperm disperse from the tissue.

5. Cryopreserving sperm: the sperm suspension was divided into 10 aliquots of 10 μ L after 3 min, Take 10 straws, carefully aspirate about 100 μ L of HTF, two 15mm of air column are spaced between 10 μ L of sperm suspension , and seal the open end of the sperm straws. Place 10 sealed straws into the freezing canister, and float them on top of the LN₂ in the cooling chamber for 10 min. After 10 min , immerse the freezing container filled with liquid N₂.

Preparation and loading of the straw for sperm cryopreservation

Cautions:

1. The cryoprotective agent equilibrated for at least 30 min at 37 $^{\circ}$ C before use, the reagent cannot be placed in -20 $^{\circ}$ C, it will cause into ingredients separate out and cannot be used.
2. When releasing the sperm, place the g-CPA drop on a hot plate at 37 $^{\circ}$ C for 3 min and swirl it in the same direction every 1 min so that the released sperm is evenly distributed in the drop.
3. When sealing sperm straws,we should be confirmed that the straws was sealed to avoid the spill of sperm , If the sperm is spilled during resuscitation from the straws, Water bath heating is prohibited and place it directly in the incubator for 10 minutes.
4. The cryoprotective agent is suitable for most strains, individual strains may have low recovery efficiency.

5. The c-TYH, HTF reagents has been required when the customer wants to transfer the embryos .
6. Before sperm cryopreservation, Mouse strain information pre-marked on the outer wall of sperm straws.
7. Ensure safety when using LN₂.