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Vitrification solution used for mouse embryo thawing (Mouse-specific)

Product	Item No.	specification	Application
0.75Su	72025	1mL	thawed embryos of mice
0.25Su	72026	1mL	thawed embryos of mice

Composition

Sucrose, Dulbecco' s Phosphate-Buffer Saline(D-PBS).

Uses and Principles

This reagent is mainly used for warming embryos of mouse embryo cryopreservation, applicable to all periods from oocyte to blastocyst, generally need to be used with KOSM or M16.

Storage conditions and expiration date

1. Storage Temperature: 4°C,Protect from light.
2. Expiration date: 6 months. Discard unused reagents after opening.

Methods

1. 0.25Su, 0.75Su and M16 were pre-warmed in incubator for more than one and a half hours for embryo thawing, and M16 dishes covered with mineral oil were prepared in advance.
2. Maintain the cryovials at room temperature for 30 seconds.
3. Pipette 900μL of 0.75Su into the cryovials, and gently blow the bottom of the cryotubel with a pipette to observe whether the liquid at the bottom is thawed. After the liquid is completely thawed, transfer the sucrose solution in the cryovials to a 35 mm dish, and then pipette 400-500μL of 0.75Su into the cryovials repeat the above operation.
4. Remove 0.25Su from the incubator and make three 50uL drops on a 35mm dish.
5. Place the dish under the microscope, collect the embryos in the 0.75Su drop, and transfer them into the first drop of 0.25Su for 1min, wash them through the other two drops of 0.25Su every 1 min. During this process,you should see that the thawed

oocytes/embryos slowly swell up due to rehydration.

6. 50 μ L M16 of three drops was pipetted into a dish, embryos were collected and washed three times in M16, and finally the number of surviving embryos was counted and placed in the mineral oil covered M16.

Cautions:

1. When embryo cryopreservation, do not cover the cryovials cap too tightly to avoid wasting time when the cap cannot be unscrewed during embryo thawing.
2. Unscrew the cryovials and immediately turn over the mouth of the tube and drain the liquid nitrogen that may have immersed in the tube.
3. Be as gentle as possible when blowing the liquid with the pipette which prevent to appear large amounts of foam.